

User Guidelines



AnyWare™

Many of our Coordinated Diagnostics Platform solutions, including Physician Portal, Advanced EMR Connectivity, Lab to Lab Connectivity, Utilization Optimization and Actionable Intelligence include ATLAS AnyWare™ technology, for sending alerts and notifications to stakeholders using a multitude of mechanisms.

Once you have one of these solutions implemented, there is nothing else that needs to be installed. You can just start setting up your desired alerts. For instance, to establish an AnyWare notification subscription in Physician Portal, follow these steps:

1. Select Notification Subscriptions from the User menu. The Notification Subscriptions page will appear.
2. Search and select the desired user to whom the notification is to be sent.
3. Search and select the desired site for which the user is to be notified. If the notification trigger (i.e., Issue Type) occurs in the selected site, the notification will be sent to the user.
4. Search and select the desired Issue Type from the lookup. This determines what event will trigger the issue notification.
5. Select the desired notification options from the Subscriptions grid.
6. Once all desired subscriptions have been added, click the Save Subscriptions button.